


Compact NMR Spectroscopy


Sub ppm resolution and high sensitivity with a benchtop system



www.act-aachen.com

www.magritek.com

- 1 Tesla magnet
- resolution < 0.04 ppm
- SNR > 30.000 per scan
- 5 mm OD sample tubes
- Clear bore aperture
- No stray magnetic field
- Fully digital electronic
- Compact Footprint
- Temperature stabilized magnet
- Shim coils up to second order



For further information contact Dr. John Trail: johntrail@magritek.com

Note - this information on this page is provided as a basis for discussion and may not reflect any final product.